







# Validation of an Indirect Fluorescent Antibody Test for Sarcocystis neurona infection in California sea lions

<u>Amalie Luneng Solli<sup>1</sup>, Cara Field<sup>2</sup>, Pádraig Duignan<sup>2</sup>, Andrea Packham<sup>1</sup>, Magdalena Plancarte<sup>1</sup>, </u> Karen Shapiro<sup>1</sup>, Devinn Sinnott<sup>1</sup>, Woutrina A. Smith<sup>1</sup> <sup>1</sup>School of Veterinary Medicine University of California, Davis <sup>2</sup>The Marine Mammal Center

### **Progress to date**

- 20 sea lion cases and 31 sea lion controls have been identified. The target sample size is 40 cases and 80 control animals.
- **Cases** are animals with suspected myopathy confirmed by necropsy and evidence of S. neurona cysts on histopathology; and will be confirmed with molecular identification.
- **Controls** are animals with a non-protozoal cause of death and no evidence on histopathology of S. neurona infection. The Kappa statistic on comparing test performance is **0.92**.

### **IFAT** titer to call a sample positive

• 1:320 titer maximizes sensitivity and specificity at 95 and 96.8%, respectively, with area under the receiver-operating curve of 0.98.

### Key Findings

- threshold for calling a sample positive.

### Discussion

- other stranded marine mammal facilities.
- environmental conditions

### Limitations

- Financial restraints in studying wildlife diseases

### **Future Directions**

- and time until symptom onset in CSL

I would like to express immense gratitude for the continued collaboration between The Marine Mammal Center and One Health Institute in addition to Lian Hortensius for providing a patient database. Financial support was provided by Students Training in Advanced Research (STAR) Program through a UC Davis School of Veterinary Medicine Endowment Fund. Image Credit © The Marine Mammal Center, Smith & Shapiro Lab.





## Discussion

 Kappa statistic compares agreement between histopathology and IFAT in classifying a case or control as positive or negative. A value of 0.92 suggests that there is good agreement between the two.

Preliminary data suggest that an IFAT titer of 1:320 is likely an appropriate

• Validation of the low-cost, ante-mortem IFAT diagnostic tool is important due to the increase of sarcocystosis cases in stranded CSL at TMMC and

Postulated factors contributing to increased *S. neurona* infections in CSL: Increase in parasite prevalence associated with changing

Shift in prey consumption based on altered historic prey availability

• Serologic antibody titer does not equal active disease • Time constraints on TMMC and UCD labs with high-volume case loads

• Risk factors associated with clinical disease and outcome • Environmental source of infection (e.g. food source, microplastics, etc.)

# Acknowledgements

The Marine Mammal Center operates under NOAA permit # 24359

Miller, M.A. et al. 2002. Evaluation of an Indirect Fluorescent Antibody Test (IFAT) for demonstration of antibodies to Toxoplasma gondii in the Sea otter (Enhydra lutris). J. of Parasitology, 88(3):549-599.

Seguel et al. 2019. Polyphasic Rhabdomyositis in California Sea Lions (Zalophus Californianus): Pathology and Potential Causes. Veterinary Pathology. 56(4):619-629.

Taylor, A. et al. 2023. Epidemiology and risk factors for Sarcocystis neurona and Toxoplasma gondii exposure in California sea lions (Zalophus californianus) from 2010-2019.