

## Clinical Diagnostic Laboratory Canine Thyroid Registry Panel Submission Form (Confidential)

General instructions: Collect blood in a plain glass or plastic tube (no anticoagulants). Place the sample in the refrigerator for 60-90 minutes to allow clotting, then centrifuge. Transfer at <u>least 2 mL of serum</u> to a plain plastic tube, labeled with owner's name and animal ID. If the specimen must be stored more than 12 hours prior to shipping, frozen storage is recommended. Ship samples FedEx overnight with cold pack in a Styrofoam box with ice packs, Monday-Thursday only. Do not send via USPS (United States Postal Service). Samples arriving unchilled or room temperature (72°F) received within 48 hours from collection time are acceptable. Samples received more than 48 hours after collection must be stored chilled or frozen and arrive at the lab at room temperature or lower. ● Non serum samples cannot be accepted. ● Hemolyzed or lipemic serum is not recommended. ● Test results will be sent by fax or mail only to the submitting veterinarian and the OFA. ● Results are not available by telephone due to confidentiality.

Enclose this request form and an OFA application form (www.offa.org) along with 2 checks: one payable to "OFA" for inclusion in the Thyroid Database and the other payable to "VMTH" to cover the cost of laboratory testing. Submissions without a check payable to the VMTH will be charged to the referring clinic. Payment to OFA is required for processing.

<u>Mailing address:</u> UC Davis VMTH, Central Laboratory Receiving, Room 1033, 1 Garrod Drive, Davis, CA 95616 <u>Questions?</u> Contact Central Lab Receiving at 530-752-VMTH (8684) or UCDVetClinicalLabs@ucdavis.edu

www.vetmed.ucdavis.edu/clinical-laboratory

\*\* Fee discounts apply to laboratory testing if more than 5 samples are submitted\*\*

Referring Clinic / Laboratory Information						
Clinic/Laboratory name						
Clinic Address						
Referring Clinician		Clinic	Clinic VMTH ID			
Phone		Fax	Fax			
	Pati	ent Inforn	nation			
Owner's last name		Regis	Registration #			
Animal name/ID		•	Breed			
Sex	Age/DOB			On therapy?		
Collection date/time		Shipr	Shipment date/time			
Sample type		# tub	# tubes submitted			
Brief History		l.				

LAB USE ONLY LAB USE ONLY

Test	Result	Reference range
Canine TSH		0.00 - 0.60 ng/mL
Thyroglobulin Autoantibody (TGAA)		Negative
Canine Free T4		0.57 – 2.93 ng/dL
# tubes: Date Temperature of sample upon receipt (°F) Type and condition of sample upon receipt	,	