


Shipping Instructions



ALL SAMPLES MUST BE SUBMITTED BY A VETERINARIAN

1. Download our submission form.
2. Fill out all Patient, Owner, Referring Veterinarian and Clinic information. Also include any pertinent past and current history of the patient. Describe what you are looking for, including clinical diagnosis.
3. Package your formalin-fixed specimen in a leak-proof container. Seal jars to ensure no liquid can leak out. Wrap containers with bubble wrap or packing peanuts and use chemical spill pads, if available.
4. Specimens should be fixed in 10% neutral buffered formalin at a 10:1 ratio. Do not ship large quantities of formalin, i.e. less than 1 liter is okay.
5. Call for current pricing 530-752-1368.

Mailing Address:

**VMTH Anatomic Pathology
VM3A, Room 1346
1 Garrod Drive
University of California, Davis
Davis, CA 95616**

**530-752-1368
530-752-7242 (FAX)**